




Evita® & More



Ventilation Pathway

The Dräger Evita® Family accompanies you along the respiration pathway and delivers high level therapy and decision support – for your daily challenges today and tomorrow.



PREVENT

- Non-invasive ventilation with sensitive inspiratory triggering and expiratory cycling-off to improve patient-ventilator synchrony

STABILIZE

- Low Flow analysis of static lung characteristics to optimize the ventilation settings
- QuickSet® and PressureLink allow to facilitate virtually any recruitment procedure
- AutoRelease in PC-APRV stabilizes the end-expiratory lung volume

WEAN

- SmartCare®/PS as an automated, knowledge-based weaning software reduces length of stay in the ICU*
- Configurable limits in SmartCare®/PS allow individual adaptations to patient situation
- Comprehensive monitoring function to observe the progress of weaning
- Variable Pressure Support to increase the natural tidal volume variability

RECOVER

- O₂-therapy supports recovery after extubation and non-invasive ventilation

*F. Lellouche et al.; Am J respir Care Med Vol 174, pp 894-900, 2006. Results are based on a European Multicenter Randomized Trial with 144 patients demonstrating improved respiratory condition, with stable hemodynamic and neurologic status, and no ARDS prior to initiating weaning

- A Dräger Evita® Infinity® V500*
- B Dräger Evita® V300
- c Dräger Evita® XL

Compatibility tested accessories and consumables that streamline your workflow and simplify your daily routine are an important part of the Evita® family product portfolio. By providing complete system solutions, Dräger focuses on increasing patient safety and reducing cost of ownership.

Please contact your local Dräger sales representative for further information. Not all devices and accessories are available worldwide.

*Infinity® Acute Care System Workstation Critical Care.



CONTENT

APPLICATIONS FOR RETROFITTING	06	BABYFLOW® PLUS	26
ACCESSORIES FOR BASIC UNIT	08	VENTILATION ACCESSORIES	29
1. Accessories for basic device	08	1. Disposable breathing circuits	29
2. Attachment parts	10	2. Reusable breathing circuits	36
3. Further accessories	10	3. Breathing system filters and HMEs	46
4. Replacement parts for application CO ₂	15	4. Catheter mounts and airway connectors	50
5. Consumables for basic device	16	5. Breathing masks	53
6. Central gas supply	19	6. Pack2Go®/Set2go	62
		7. Resuscitators	63
BREATHING GAS HUMIDIFIERS	20		
1. Humidification with F&P MR 850	20		
MEDICATION NEBULIZER	24		
1. Pneumatic medication nebulizer	24		
2. Aerogen Aeroneb Pro	25		

APPLICATIONS FOR RETROFITTING

APPLICATIONS FOR RETROFITTING

Retrofit Kits / Applications

(when ordering state: country-specific features)

8416200

Hardware Applications

Nebulizer

Nurse Call

CO₂ monitoringAssembly plate 38 mm (1.5 inch) pole
for mounting of ventilation unit

Trolley 2 – 90 cm

Supporting frame with DIN rail

Supporting frame with DIN rail and bail

Cockpit holder

A Gas supply unit GS500

Power supply unit PS500

B Transport Supply Unit

Transport coupling for Transport Supply Unit

Gas cylinder holder

C Bed coupling

Bed side rail clamp for bed coupling

Software Applications

D Airway Pressure Release Ventilation

Low Flow PV Loop

AutoFlow/Volume Guarantee

Mandatory Minute Ventilation

Automatic Tube Compensation

Proportional Pressure Support

Monitoring Plus (advanced trend parameter, tabular trends,
graphical trends, data export and advanced loops)

Neonatal Ventilation

E SmartCare®/PS

F Non-Invasive Ventilation

G Smart Pulmonary View

Pressure Support Ventilation

Variable Pressure Support

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A

D-11094-2011




B

D-43489-2012




C

D-43526-2012




D

D-373-2017




E

D-374-2017



F

D-375-2017



G

D-376-2017



APPLICATIONS FOR RETROFITTING

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

Evita® XL kits for retrofitting	
Order number for the kits listed below	8414240
CapnoPlus kit	
Nurse call kit	
Evita® Remote kit	
Communication / Evita®Link kit	
NeoFlow kit	
Mask ventilation (NIV) kit	
Evita® XL 5.10 Evita® XL software kit	
6.0 Evita® XL software kit	
A SmartCare® kit; not for USA	
SmartCare® Update 1.1	
SmartCare® Capno package kit	
Lung protection package kit	
O ₂ -therapy kit	
Proportional Pressure Support (PPS) kit; not for USA	
Software 7.06, NIV Plus	
Plug for nurse call connection	1846248
Accessories for nurse call option.	
B CapnoSmart CO₂ main flow sensor	6871500
Replacement parts for SmartCare® option (not for USA).	
For patient application see O2Star Portfolio.	

A

MT-0619-2008



B

MT-0686-2007



ACCESSORIES FOR BASIC UNIT

ACCESSORIES FOR BASIC UNIT

1. ACCESSORIES FOR BASIC DEVICE

ACCESSORIES FOR TROLRIES EVITA® INFINITY® V500 AND EVITA® V300

A	Gas cylinder holder	G93110
B	Bedside rail clamp for bed coupling	8419801
C	Rail clamp to mount F&P humidifier	8411074

ACCESSORIES FOR MEDICAL COCKPIT NOT MOUNTED ON TROLLEY

	Module M26 for mounting of Infinity® Medical Cockpit (when ordering state: country-specific features)	G96500
	System cable, length 0.70 m (27.6 inch)	MS18782
	System cable, length 2.80 m (9.2 ft)	MS18781

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



ACCESSORIES FOR BASIC UNIT

TROLLEY AND ACCESSORIES FOR TROLLEY EVITA® XL

A	EvitaMobil trolley (high); not for USA	8411950
	EvitaMobil trolley (low); not for USA	8411965
	EvitaMobil with extension (USA)	8414455

ACCESSORIES FOR TROLRIES EVITA® XL

B	Bottle holder set, EvitaMobil	8411970
C	DC battery 12 V/17 Ah (box of 1; 2 required)	1843303
	Battery cable kit for EvitaMobil	8411822
D	Pressure reducer OxyLine, customized configuration of pressure reducers	MP00100
	OxyLine's modular design means you can accurately choose the pressure reducer you need from the following components:	
	– O ₂ or AIR pressure reducer with 200 or 300 bar	
	– Manometer, vertical or horizontal	
	– Flow regulator settings (1, 2, 3, 4, 6, 8, 12 and 15 L/min) and 9/16" sterile water connector	
	– Flowmeter for up to 16 L/min with connector for Dräger reusable humidifier, sterile water connector and 6 mm (0.2 in) nozzle for an oxygen mask	
	– Self-closing NIST connector	
	– Centralized gas supply coupling, e.g. to connect up a mains-operated ventilation device	
	– Ejector for short-time suction of eg. secret	
E	O₂-NIST/NIST Connecting hose 1.2 m	M36053
	AIR-NIST/NIST Connecting hose 1.2 m	M36055
	A choice of additional hoses for connecting the pressure reducer with the Evita XL: 8601697.	

COMPRESSOR EVITA® XL

F	Breathing air compressor for supply with compressed air, can be integrated in the EvitaMobil	8414350
	Standalone base for compressor	8411520
	Air hose, 0.8 m, NIST neutral	8411517
	Air hose, 0.8 m, NIST ISO 32	8411538
	Air hose, 5 m, NIST ISO 32	8411541
	Air hose, 5 m, NIST neutral	8411519
	Air hose, 0.8 m, DISS	8414454

Not all articles are available worldwide.

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



ACCESSORIES FOR BASIC UNIT

2. ATTACHMENT PARTS

A	Instrument tray, complete (EvitaXL only)	8414828
	Side holder set	8411615

STORAGE

B	CS connection hose holder	2M85446
C	Hook	M26349
	Hook for accessories, e.g. breathing bag.	

3. FURTHER ACCESSORIES

	Coarse dust filter for ventilator, for patient part heating, as of 2005	8415978
	Cover for flow sensor, blue, for patient part heating, up to 2005	8414714
	Cover for flow sensor, Plexiglas, for patient part heating, up to 2005	MP01003
	Flow sensor flap; for patient part heating, as of 2005	8415973
	ILV cable, Evita® 4/XL	8411794
	ILV cable, Evita® 1/2	8411793
	Monitor bracket kit, Evita®, with trolley	8415787
	Monitor bracket kit, Evita®, without trolley	8415850
	Printer cable	8306489

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



ACCESSORIES FOR BASIC UNIT


MEDIBUS cable	8416326	■	■	■
Software Protocol, MEDIBUS for V and VN ventilators	9039527	■	■	
Software Protocol, MEDIBUS.X for V and VN ventilators	9052608	■	■	
Protocol Definition, Dräger RS 232 MEDIBUS	9028258	■	■	
Instructions for Use, MEDIBUS for Intensive Care Devices	9028329			■
Instructions for Use, RS 232 MEDIBUS Protocol Definition	9028320			■
USB flash drive	8416347	■	■	

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

ACCESSORIES FOR BASIC UNIT

A	Test lung SelfTestLung, latex-free, reusable, max. 1 L (33.81 fl. oz.)	MP02400
	Spare parts set SelfTestLung, with locking brace and attachment bolt	MP02401
B	Test lung NeoTestLung, reusable, N, latex-free, bellows	8409742
C	Demo lung to simulate neonatal lungs, with suitcase	8418077

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL






ACCESSORIES FOR BASIC UNIT

HINGED ARMS

Keeping everything close at hand. To leave your hands free to care for the patient.

A	Hinged arm Quickstop 1, Length 106.7 cm, max. load 2 kg	2M85705			
B	Hinged arm Quickstop 2, w/ universal ball holder system for optional accessories	2M85706			
C	Hinged arm 240°, 3 hinges, fixable, ball holder system, with hose clamps	8409609			
D	Hinged arm Infinity® ACS, for positioning and fixation of breathing hoses	MP00690			

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



ACCESSORIES FOR BASIC UNIT



A	Hose clip, adult, neo., for ball mounting, for standard breathing hoses	8409841
B	Hose clip double, for standard breathing hoses, adult, neo., 2 pcs.	8406138
C	Holder for Y-piece, with ball mounting	8409746
	Hose clip single	8408356
D	Universal hose clip, for hinged arm Infinity® ACS (MP00690)	MP00692
E	Hose clip double, gray, for hinged arm Infinity® ACS (MP00690)	MP00693
F	CS connection hose holder	2M85446

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

ACCESSORIES FOR BASIC UNIT

4. REPLACEMENT PARTS FOR APPLICATION CO₂

A	CO ₂ mainstream sensor	6871950
B	CapnoSmart CO ₂ main flow sensor Replacement parts for SmartCare® option (not for USA).	6871500
CO ₂ CUVETTES, DISPOSABLE		
C	① CO ₂ cuvette disposable, pediatric, 10 pcs.	MP01063
	② CO ₂ cuvette disposable, adult, 10 pcs.	MP01062
CO ₂ CUVETTES, REUSABLE		
D	① CO ₂ cuvette reusable, pediatric, for CO ₂ sensors 6871500, 6870300, 6871950	6870280
	② CO ₂ cuvette reusable, adult, for CO ₂ sensors 6871500, 6870300, 6871950	6870279
ACCESSORIES		
	Calibration set, with test connection and test gas cylinder	8412710
	Test gas cylinder 5 Vol.% CO ₂ , 95 Vol. % N ₂	6850435
E	Parking holder for CO ₂ sensor	8416352
	Parking holder for CO ₂ sensor	8412840

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

ACCESSORIES FOR BASIC UNIT

5. CONSUMABLES FOR BASIC DEVICE

A	Expiratory valve for Evita® 2 dura, Evita® 4 edition, Evita® XL, disposable, 10 pcs.	8414776
	Expiratory valve for Evita® 2 dura, Evita® 4 edition, Evita® XL, reusable	8410580
B	Expiratory valve for Evita® 2 dura, Evita® 4 edition, Evita® XL, disposable, 10 pcs.	MP02600
	Membrane for Evita® 2 dura, Evita® 4 edition, Evita® XL, for expiratory valve	8412015
C	Water trap for Evita® 2 dura, Evita® 4 edition, Evita® XL, reusable	8413125
D	Expiratory valve for Evita® Infinity® V500, Evita® V300, disposable, Infinity® ID, 10 pcs.	MP01060
	Expiratory valve for Evita® Infinity® V500, Evita® V300, reusable	8416750
E	Infinity® ID Expiratory filter, disposable, 20 pcs.	MP01780
F	Expiratory filter, disposable, 20 pcs.	MP01781
COMPATIBLE WITH THE EXPIRATORY VALVE FOR ADULTS/PEDIATRIC PATIENTS		
	Replacement diaphragm for reusable expiratory valve	8413661
	Infinity® ID Inspiratory valve, reusable	8415195
G	Infinity® ID Expiratory valve, reusable, neonates/pediatric patients	8415270
	Ambient air filter/Filter for fan, set of 10 pcs.	8416349
	Ambient air filter/Filter for fan, set of 10 pcs.	8412384

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

ACCESSORIES FOR BASIC UNIT

FLOW SENSORS BY DRÄGER

Hot-wire anemometry: Accurate and intelligent flow measurement


Spirolog®, neonatal flow sensors and SpiroLife® are based on proven hot-wire anemometry techniques. Due to extremely fast and accurate measurement results, as well as extremely low inspiratory and expiratory resistances, it is clearly superior to other methods. You only need one flow sensor for all patients. Neonatal flow sensors offer customized solutions for the sole exception of the specific needs of neonates in intensive care. The sterilizable neonatal flow sensors and SpiroLife® combine the advantages of hot-wire technology with the possibility of sterilizability. Therefore, they integrate seamlessly into your hospital's sterilization processes.

A	Flow sensor Spirolog®, PC/ABS, disinfectable, 5 pcs.	8403735	■
B	Flow sensor SpiroLife®, PPSU, disinfectable and sterilizable	MK01900	■
C	Flow sensor Infinity® ID, PC/ABS, disinfectable, 5 pcs.	6871980	■ ■
D	Neonatal flow sensor with Y-piece, PSU, disinfectable and sterilizable	8410185	■ ■
E	Neonatal flow sensor ISO 15, PSU, disinfectable and sterilizable	8411130	■ ■ ■
F	Insert for neonatal flow sensor, PSU, disinfectable and sterilizable, 5 pcs.	8410179	■ ■ ■
	Flow sensor cover	8417163	■

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A

D-21179-2015



B

D-7360-2011



C

D-7483-2016



D

MT-351-2003



E

D-9121-2009



F

MT-0510-2007



ACCESSORIES FOR BASIC UNIT

OXYGEN SENSOR BY DRÄGER

Oxygen for Life

Oxygen sensors have the important task of measuring and monitoring the inspiratory and/or expiratory O₂ concentration in anesthesia and respiratory devices. Oxygen sensors from Dräger satisfy the highest national and international requirements. Furthermore, they have been extensively tested and fine-tuned for proper functional performance with Dräger devices and systems. Dräger's years of experience in the field of sensors, and their high quality has even convinced NASA. Dräger O₂ sensors are electrochemical transducers used to measure O₂ partial pressure. They operate on the principle of a galvanic cell.

Examples of locations where the sensors are used include the Fabius® family, Evita® family and the Primus®.

To take into account the high value of oxygen measurement, we only recommend using O₂ sensors from Dräger.

A	O₂ sensor capsule, for mainstream measurement	6850645
----------	---	----------------

CABLES

B	Flow sensor cable, for neonatal flow sensor, e.g., Babylog	8409626
C	Flow sensor cable for V500 with Neo option and VN500	8416600

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A

MT-423-2002



B

MT-2672-2004



C

D-9315-2011



ACCESSORIES FOR BASIC UNIT

6. CENTRAL GAS SUPPLY

A	Dräger supply hoses, selectable by the adjacent configuration number	8601697
	Design, color coding and connection thread of all Dräger supply hoses comply with the respective legal regulations in the indicated countries.	

FILTER FOR CENTRAL SUPPLY CONNECTING HOSES

High Pressure Line Filter for medical air	MP03901
High Pressure Line Filter for oxygen	MP03902
Filter element with O-ring	MP03903
Filter bowl with O-ring	MP03904

ADAPTER HOSES

Adapter O₂ NIST for Australia	M34446
Adapter AIR NIST for Australia	M34448

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A

D-2286-2014




BREATHING GAS HUMIDIFIERS

BREATHING GAS HUMIDIFIERS

1. HUMIDIFICATION WITH F&P MR 850

	Humidifier basic unit F&P MR 850 (when ordering state: country-specific features)	8414144			
A	Temperature probe/flow sensor for F&P MR 850, for all breathing circuits, 1.1 m (3.6 ft)	8414989			
A	Temperature probe/flow sensor for F&P MR 850, for all breathing circuits, 1.5 m (4.92 ft)	8414966			
B	Hose heater adapter for F&P MR 850, for reusable breathing circuit	8414968			
C	Hose heater adapter for F&P MR 850, for disposable breathing circuit	8414992			
D	Hose heater adapter for F&P MR 850, for double-heated disposable breathing circuit	MP01006			
E	Draw wire for hose heater 1.5 m, for reusable breathing circuit	8411050			


EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



BREATHING GAS HUMIDIFIERS

REQUIRED ACCESSORIES FOR F&P MR 850



A	Humidifier chamber for F&P MR 290, disposable, 10 pcs.	8418282	■	■	■
B	Humidifier chamber F&P MR 340, reusable, N/P	8411047	■	■	■
C	Humidifier chamber F&P MR 370, reusable, A	8412217	■	■	■
D	Humidifier chamber VentStar®, disposable, 10 pcs.	MP00234	■	■	

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

BREATHING GAS HUMIDIFIERS

HOLDER

A	Rail clamp to mount F&P humidifier	8411074
B	Swivel-mounted humidifier holder for trolley	G93111
C	Holder for humidifier, lateral standard rail	8416325
D	Universal holder with standard rail	G93140
E	Holder for F&P humidifier U-shaped, for trolley, with rail clamp	MP01002
	Humidifier holder	8411956

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

BREATHING GAS HUMIDIFIERS

HOSES

A	Hose heater, A, 1.1 m	8411045
B	Hose heater for F&P MR 700/730 and MR 850, for reusable breathing circuit, 1.1 m	8414993

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A

MT-2264-2004



B

MT-2266-2004



MEDICATION NEBULIZER

MEDICATION NEBULIZER
1. PNEUMATIC MEDICATION NEBULIZER

A	Medication nebulizer for Evita® 2 dura, Evita® 4 edition, Evita® XL	8412935
B	Connector set for medication nebulizer, for neonatal breathing circuit	8411031

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



MEDICATION NEBULIZER

2. AERONEB® PRO

Aeroneb Pro works with an electronic micropump. The nebulizer can be used with Evita® 2 dura, Evita® 4 edition and Evita® XL.

Aeroneb Pro is ideally suited to nebulizer therapy requirements for inline mechanical ventilation and it accommodates all patient age ranges from neonates to adults. It operates without changing patient ventilator parameters, without heating or degrading medications and it can be refilled without interrupting ventilation.

A	Medication nebulizer Aeroneb Pro	MP01010			
B	Nebulizer unit with filler cap for Aeroneb Pro	MP01030			
C	T-piece with plug for Aeroneb Pro, A, ISO 22 mm (0.87 in), 5 pcs.	MP01031			
D	T-piece set with plug for Aeroneb Pro, P, 15 mm (0.59 in), w/ adapter 10 mm, 5 pcs.	MP01032			
E	T-connector for Aeroneb Pro, 12 mm (0.47 in), with plug, 5 pcs.	MP01037			
F	Controller cable for Aeroneb Pro	MP01036			
G	Filler cap for Aeroneb Pro, silicone, 5 pcs.	MP01035			

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



BABYFLOW® PLUS


BABYFLOW® PLUS

BabyFlow® plus – The superior non-invasive respiratory support system, designed to maximise performance and patient comfort.

BABYFLOW® PLUS

A	BabyFlow® plus, disposable, size S, 20 pcs.	MP03700	■
	BabyFlow® plus, disposable, size M, 20 pcs.	MP03701	■
	BabyFlow® plus, disposable, size L, 20 pcs.	MP03702	■
B	Demo package BabyFlow® plus, disposable, size S	MP03703	■
	Demo package BabyFlow® plus, disposable, size M	MP03704	■
	Demo package BabyFlow® plus, disposable, size L	MP03705	■
C	Soft connektor N, disposable, 20 pcs.	MP03826	■

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



BABYFLOW® PLUS

STANDARD PRONGS

Soft anatomically formed nasale silicone prongs, disposable.

A	① Prong BabyFlow®, disposable, size XS, 10 pcs.	8418415	■	■	■
	② Prong BabyFlow®, disposable, size S, 10 pcs.	8418605	■	■	■
	③ Prong BabyFlow®, disposable, size M, 10 pcs.	8418416	■	■	■
	④ Prong BabyFlow®, disposable, size L, 10 pcs.	8418531	■	■	■
	⑤ Prong BabyFlow®, disposable, size XL, 10 pcs.	8418417	■	■	■
B	Prong BabyFlow®, disposable, size XXL, 10 pcs.	8418617	■	■	■

WIDE PRONGS




Prong BabyFlow®, disposable, size 3, 10 pcs.	8418603	■	■	■
Prong BabyFlow®, disposable, size 4, 10 pcs.	8418607	■	■	■
Prong BabyFlow®, disposable, size 5, 10 pcs.	8418609	■	■	■
Prong BabyFlow®, disposable, size 6, 10 pcs.	8418611	■	■	■

MASKS

Nasal masks of silicone and TPE complete the portfolio and increase the comfort for the patient!

C	① Neo-mask, silicone, disposable, size S, 10 pcs.	8418491	■	■	■
	② Neo-mask, silicone, disposable, size M, 10 pcs.	8418490	■	■	■
	③ Neo-mask, silicone, disposable, size L, 10 pcs.	8418619	■	■	■
D	Neo-mask, TPE, disposable, size 4, 10 pcs.	MP01494	■	■	■
E	Neo-mask, TPE, disposable, size 5, 10 pcs.	MP01495	■	■	■
F	Neo-mask, TPE, disposable, size 6, 10 pcs.	MP01496	■	■	■

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



BABYFLOW® PLUS

HEADGEARS

The new headgear enable the customer to choose the most comfortable fit and combination for each of their tiny patients.

- A** Headgear, disposable, size XXS,
17 – 20 cm (6.7 – 7.9 in) head circumference, 5 pcs. MP03880

Headgear, disposable, size XS,
20 – 23 cm (7.9 – 9.1 in) head circumference, 5 pcs. MP03881

- B** Headgear, disposable, size S,
23 – 27 cm (9.1 – 10.6 in) head circumference, 5 pcs. MP03882

Headgear, disposable, size M,
27 – 31 cm (10.6 – 12.2 in) head circumference, 5 pcs. MP03883

Headgear, disposable, size L,
31 – 37 cm (12.2 – 14.6 in) head circumference, 5 pcs. MP03884

Headgear, disposable, size XL,
37 – 42 cm (14.6 – 16.5 in) head circumference, 5 pcs. MP03885

CAPS

The disposable and breathable microfiber caps (95% cotton, 5% spandex) are available in seven sizes (color-coded)

- C** **①** Cap BabyFlow®, disposable, size XS, white,
18 – 20.5 cm (7.1 – 8.1 in) head circumference, 5 pcs. 8418533

Cap BabyFlow®, disposable, size S, yellow,
20.5 – 23 cm (8.1 – 9.1 in) head circumference, 5 pcs. 8418534

Cap BabyFlow®, disposable, size M, red,
23 – 25.5 cm (9.1 – 10 in) head circumference, 5 pcs. 8418535

- ②** Cap BabyFlow®, disposable, size L, green,
25.5 – 28 cm (10 – 11 in) head circumference, 5 pcs. 8418536

- ③** Cap BabyFlow®, disposable, size XL, orange,
28 – 30.5 cm (11 – 12 in) head circumference, 5 pcs. 8418537

Cap BabyFlow®, disposable, size XXL, blue,
30.5 – 33 cm (12 – 13 in) head circumference, 5 pcs. 8418538

Cap BabyFlow®, disposable, size XXL+, white/black,
33 – 36 cm (13 – 14.2 in) head circumference, 5 pcs. 8418539

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



VENTILATION ACCESSORIES


VENTILATION ACCESSORIES

1. DISPOSABLE BREATHING CIRCUITS

INFINITY® ID PRODUCTS FOR ADULTS

A	Breathing circuit Infinity® ID, latex-free, disposable, basic, 1.8 m (5.91 ft), 25 pcs.	MP01350
B	Breathing circuit Infinity® ID, latex-free, disposable, basic, 2.5 m (8.2 ft), 10 pcs.	MP01348
C	Breathing circuit Infinity® ID, latex-free, disposable, water trap, 1.8 m (5.91 ft), 25 pcs.	MP01326
D	Breathing circuit Infinity® ID, latex-free, disposable, coaxial, 1.6 m (5.25 ft), 20 pcs.	MP00342
E	Breathing circuit Infinity® ID, latex-free, disposable, smoothbore, 1.6 m (5.25 ft), 25 pcs.	MP00340
F	Breathing circuit Infinity® ID, latex-free, disposable, smoothbore, 3 m (9.8 ft), 10 pcs.	MP01325
G	Breathing circuit Infinity® ID, latex-free, disposable, extendable, 2.2 m (7.22 ft), 25 pcs.	MP01347

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



DISPOSABLE BREATHING CIRCUITS

A Breathing circuit Infinity® ID, latex free, disposable, coaxial, 1.9 m (6.2 ft), 20 pcs. MP03375

B Breathing circuit Infinity® ID, latex free, disposable, coaxial, 2.3 m (7.55 ft), 20 pcs. MP03377

C Breathing circuit Infinity® ID, latex free, disposable, coaxial, with water traps, 2.1 m (6.89 ft), 10 pcs. MP00378

INFINITY® ID PRODUCTS FOR PEDIATRICS

D Breathing circuit Infinity® ID, latex-free, disposable, P, basic, 1.8 m (5.91 ft), 25 pcs. MP01340

E Breathing circuit Infinity® ID, latex-free, disposable, P, water trap, 1.8 m (5.91 ft), 25 pcs. MP01341

INFINITY® ID PRODUCTS FOR NEONATES

F Breathing circuit Infinity® ID, latex-free, disposable, N, heated, 1.2 m (3.9 ft), 10 pcs. MP00343




EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



DISPOSABLE BREATHING CIRCUITS



PRODUCTS WITHOUT INFINITY® ID FOR ADULTS

A	Breathing circuit Vent Set, latex-free, disposable, basic, 1.5 m (4.92 ft), 25 pcs.	MP00300
B	Breathing circuit VentStar®, latex-free, disposable, basic, 1.8 m (5.91 ft), 25 pcs.	MP00349
C	Breathing circuit VentStar®, latex-free, disposable, basic, 2.5 m (8.2 ft), 10 pcs.	MP00350
D	Breathing circuit Vent Set, latex-free, disposable, extendable, 1.8 m (5.91 ft), 25 pcs.	MP00305
E	Breathing circuit VentStar®, latex-free, disposable, extendable, 2.2 m (7.22 ft), 25 pcs.	MP00355

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

DISPOSABLE BREATHING CIRCUITS


A	Breathing circuit Vent Set, latex-free, disposable, coaxial, 1.5 m (4.92 ft), 20 pcs.	MP00315
B	Breathing circuit Vent Set, latex-free, disposable, coax, with holder, 1.5 m (4.92 ft) 20 pcs.	MP00316
C	Breathing circuit VentStar®, latex-free, disposable, water trap, w/o LL, 1.8 m (5.91 ft), 25 pcs.	MP00334
D	Breathing circuit VentStar®, latex-free, disposable, water trap, 1.8 m (5.91 ft), 25 pcs.	MP00337
E	Breathing circuit VentStar®, latex-free, disposable, smoothbore, for MRI, 3 m (9.8 ft), 10 pcs.	MP00327
F	Breathing circuit VentStar®, latex-free, disposable, coaxial, 1.9 m (6.2 ft), 20 pcs.	MP03374
G	Breathing circuit VentStar®, latex-free, disposable, coaxial, with water traps, 2.1 m (6.89 ft), 10 pcs.	MP00373

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

DISPOSABLE BREATHING CIRCUITS

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A	<u>VentStar® Helix, dual heated, 10 pcs.</u>	<u>MP02606</u>
B	<u>VentStar® Helix, heated, 10 pcs.</u>	<u>MP02607</u>
C	<u>Humidifier chamber, 10 pcs.</u>	<u>MP02605</u>
	<u>Breathing circuit Fisher & Paykel RT212, disposable, inspiratory heated, 10 pcs.</u>	<u>8414991</u>





DISPOSABLE BREATHING CIRCUITS

DISPOSABLE BREATHING CIRCUITS FOR PEDIATRICS

A	Breathing circuit VentStar®, latex-free, disposable, P, basic, w/o LL, 1.5 m (4.92 ft), 25 pcs.	MP00338
B	Breathing circuit VentStar®, latex-free, disposable, P, basic, 1.8 m (5.91 ft), 25 pcs.	MP00351
C	Breathing circuit VentStar®, latex-free, disposable, P, water trap, 1.8 m (5.91 ft), 25 pcs.	MP00361
D	Breathing circuit VentStar®, latex-free, disposable, P, water trap, w/o LL, 1.8 m (5.91 ft), 25 pcs.	MP00362


EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



DISPOSABLE BREATHING CIRCUITS

DISPOSABLE BREATHING CIRCUITS
FOR NEONATES

A	Breathing circuit VentStar®, latex-free, disposable, N, heated, 1.2 m (3.9 ft), 10 pcs.	MP00308
B	Breathing circuit VentStar®, latex-free, disposable, N, heated, w/o chamber, 1.2 m (3.9 ft), 10 pcs.	MP00309
C	Breathing circuit VentStar® Helix, latex-free, disposable, N, dual heated, 1.7 m, 10 pcs.	MP02650
D	Breathing circuit VentStar® Helix, latex-free, disposable, N, heated, 1.7 m, 10 pcs.	MP02608
E	Humidifier chamber VentStar®, disposable, 10 pcs.	MP00234
	Breathing circuit Fisher & Paykel RT225, disposable, N, inspiratory heated, 10 pcs.	8414994

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL


REUSABLE BREATHING CIRCUITS

2. REUSABLE BREATHING CIRCUITS

BLUESET VENTILATION

A	BlueSet Ventilation HME (P/A)	8412860
	Consists of:	
	① Breathing hose, silicone, reusable, A, 1.5 m (4.92 ft), connector 22 mm/22 mm (0.87/0.87 in)	2166046
	② Y-piece, reusable, for temperature probe adapter 900 MR 532	8412888


EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



REUSABLE BREATHING CIRCUITS


BLUESET HEATED

A	BlueSet heated (P/A)	8414986
	Consists of:	
1	Breathing hose, silicone, reusable, A, 1.1 m (3.61 ft), connector 22 mm/22 mm	2166038
2	Hose clip, 10 pcs.	8403566
3	Breathing hose, silicone, reusable, A, 60 cm (1.97 ft), connector 22 mm/22 mm	2165627
4	Connector for breathing hose, reusable, A	8403076
5	Water trap, reusable, A	8404985
6	ISO double connector, reusable, 22 mm/22 mm (0.87/0.87 in)	M25647
7	Connector set for endotracheal tube, reusable, size 6 mm – 12.5 mm, 12 pcs.	8403685
8	Connector for endotracheal tube, reusable, size 12.5 mm	M23840
9	Cap, reusable, 11 mm	8401644
10	Corrugated hose, reusable, catheter mount	8402041
11	Y-piece, reusable, for temperature probe adapter 900 MR 532	8412888
12	Temperature probe adapter 900 MR 532, for F&P humidifiers	MP01007
13	Elbow, reusable	M25649
14	Hose heater for F&P MR 700/730 and MR 850, for reusable breathing circuit, 1.1 m (3.61 ft)	8414993
15	Breathing hose, silicone, reusable, A, 35 cm, connector 22 mm/22 mm (0.87/0.87 in)	2165619



REUSABLE BREATHING CIRCUITS


A	BlueSet heated (N)	8414987
	Consists of:	
❶	Breathing hose, silicone, reusable, P, 1.1 m (3.61 ft), connector 10 mm/22 mm	2165651
❷	Breathing hose, silicone, reusable, P, 60 cm (1.97 ft), connector 10 mm/10 mm	2165848
❸	Corrugated hose, silicone, 13 cm, Ø 10 mm (0.39 in), catheter mount, 2 pcs.	8409634
❹	Breathing hose, silicone, reusable, P, 40 cm (15.8 in), connector 10 mm/22 mm	2165856
❺	Water trap, reusable, N	8409627
❻	Hose heater for F&P MR 700/730 and MR 850, for reusable breathing circuit, 1.1 m (3.61 ft)	8414993
❼	Cap, reusable, 7 mm	8401645
❽	Y-piece, reusable, P, angled, with suction opening	8403075
❾	Connector for endotracheal tube, reusable, size 11 mm	M19351
❿	T-piece connector, reusable	8411044
⓫	Double connector, outer diameter 11 mm	8409897
⓬	Test lung NeoTestLung, reusable, N, latex-free, bellows	8409742
⓭	Breathing hose, silicone, reusable, P, 60 cm (1.97 ft), connector 10 mm/22 mm (0.39/0.87 in)	2165821



EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

REUSABLE BREATHING CIRCUITS

A	BlueSet heated Incubator (N)	8411041
	Consists of:	
1	Breathing hose, silicone, reusable, P, 1.1 m, connector 10 mm/22 mm	2165651
2	Breathing hose, silicone, reusable, P/N, 1 m, Ø 10 mm	8403080
3	Breathing hose, silicone, reusable, P/N, 35 cm, Ø 10 mm	8403070
4	Test lung NeoTestLung, reusable, N, latex-free, bellows	8409742
5	Water trap, reusable, N	8409627
6	Elbow, reusable	8418289
7	Double connector, outer diameter 11 mm	8409897
8	Connector for endotracheal tube, reusable, size 11 mm	M19351
9	T-piece connector, reusable	8411044
10	Hose heater, A, 1.1 m	8411045
11	Corrugated hose, silicone, N, 20 cm, catheter mount	8410709
12	Breathing hose, silicone, reusable, N, 50 cm, Ø 10 mm	8418458



EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

REUSABLE BREATHING CIRCUITS

Silicone breathing hose (adult) with
standard sleeves 22 mm/22 mm (0.9 in/0.9 in)

A	Breathing hose, silicone, reusable, 0.35 m (13.7 in), connector 22 mm/22 mm (0.87/0.87 in)	2165619
B	Breathing hose, silicone, reusable, 60 cm (1.97 ft), connector 22 mm/22 mm (0.87/0.87 in)	2165627
C	Breathing hose, silicone, reusable, 1.1 m (3.6 ft), connector 22 mm/22 mm (0.87/0.87 in)	2166038
D	Breathing hose, silicone, reusable, 1.5 m (4.9 ft), connector 22 mm/22 mm (0.87/0.87 in)	2166046

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A

D-36204-2011



B

D-36208-2011




C

D-36211-2011



D

D-36225-2011






REUSABLE BREATHING CIRCUITS

Silicone breathing hose (pediatric) with standard sleeves with one 10 mm (0.4 in) and one 22 mm (0.9 in) sleeve and 10 mm (0.4 in) inner diameter

A	Breathing hose, silicone, reusable, P, 40 cm (15.8 in), connector 10 mm/22 mm	2165856			
B	Breathing hose, silicone, reusable, P, 60 cm (1.97 ft), connector 10 mm/22 mm	2165821			
C	Breathing hose, silicone, reusable, P, 1.1 m (3.61 ft), connector 10 mm/22 mm	2165651			
D	Breathing hose, silicone, reusable, P, 1.5 m (4.92 ft), connector 10 mm/22 mm, Ø 10 mm	2165961			



EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



REUSABLE BREATHING CIRCUITS

Silicone breathing hose (pediatric) with standard sleeves 10 mm (0.4 in) on both sides and 10 mm (0.4 in) inner diameter		
A	Breathing hose, silicone, reusable, P, 60 cm (1.97 ft), connector 10 mm/10 mm	2165848
B	Breathing hose, silicone, reusable, P, 1.1 m (3.61 ft), connector 10 mm/10 mm	2165600
Silicone breathing hose (pediatric) with standard sleeves 22 mm (0.9 in) on both sides and 10 mm (0.4 in) inner diameter		
C	Breathing hose, silicone, reusable, P, 1.1 m (3.61 ft), connector 22 mm/22 mm, Ø 10 mm	2165678

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



VENTILATION ACCESSORIES

Y-PIECES AND CONNECTORS

A	Y-piece, reusable, for temperature probe adapter 900 MR 532	8412888			
B	Temperature probe adapter 900 MR 532, for F&P humidifiers	MP01007			
C	Y-piece, reusable, for Dräger temperature probe	8405435			
D	Y-piece, reusable, P, angled, with suction opening	8403075			
	Cap, reusable, 7 mm (0.28 in)	8401645			
E	Elbow, reusable	M25649			

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A

MT-1321-2008

**B**

MT-4036-2006

**C**

MT-2688-2004

**D**

MT-2726-2004

**E**

MT-2688-2004



VENTILATION ACCESSORIES

A	ISO double connector, reusable, 22 mm/22 mm (0.87/0.87 in)	M25647
	Connector for breathing hose, reusable, A	8403076
	T-piece connector, reusable	8411044
	Double connector, outer diameter 11 mm	8409897
B	Water trap	8404985
	For use between the ventilator hoses.	
C	Water trap, reusable, (N/P)	8409627


EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



VENTILATION ACCESSORIES

	Connector for endotracheal tube (catheter connector), size 12.5, set of 10 pcs.	<u>M23840</u>
A	Connector set for endotracheal tube, reusable, size 6 mm – 12.5 mm, 12 pcs.	<u>8403685</u>
	Connector for endotracheal tube, size 11	<u>M19351</u>
	Connector for tracheal tube, size 11	<u>M19490</u>
	Set of caps, set of 5 pcs.	<u>8402918</u>
	Corrugated hose	<u>8402041</u>
	Corrugated hose, flexible, length 0.14 m (5.5 in)	<u>8403333</u>
	Corrugated hose, flexible, length 0.13 m (5.3 in)	<u>8409634</u>
	Hose clip	<u>8403566</u>
	Sleeve	<u>M25647</u>
	Temperature sensor	<u>8405371</u>
B	Thermometer	<u>2M13259</u>

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



BREATHING SYSTEM FILTERS AND HMEs

3. BREATHING SYSTEM FILTERS AND HMEs

ELECTROSTATIC FILTERS AND HMEs

The TwinStar® filter/HME from Dräger is a combination of an efficient HME and a high-performance breathing system filter, thereby significantly contributing to infection prophylaxis in ventilation treatment.

- Combination of efficient HME and high-performance breathing system filter
- High retention rates for bacteria and viruses
- High return of heat and moisture
- Reduced condensation in the breathing circuit
- Minimal workload and low cost
- Sampling connector with tethered cap
- Sampling port in convenient 45° angle
- Transparent housing for visual control
- Standardized connectors for safe connection to other components
- Clear labeling and blue color coding for quick identification

A	Electrostatic filter and HME TwinStar® 90, disposable, 90 mL (3.04 fl. oz.) deadspace, 50 pcs.	MP01800	■	■	■
B	Electrostatic filter and HME TwinStar® 65A, disposable, 65 mL (2.2 fl. oz.) deadspace, angled, 50 pcs.	MP01810	■	■	■
C	Electrostatic filter and HME TwinStar® 55, disposable, 55 mL (1.86 fl. oz.) deadspace, 50 pcs.	MP01805	■	■	■
D	Electrostatic filter and HME TwinStar® 25, disposable, 25 mL (0.85 fl. oz.) deadspace, 50 pcs.	MP01815	■	■	■
E	Electrostatic filter and HME TwinStar® 10A, disposable, 10 mL (0.34 fl. oz.) deadspace, angled, 50 pcs.	MP01825	■	■	■
F	Electrostatic filter and HME TwinStar® 8, disposable, 8 mL (0.27 fl. oz.) deadspace, 50 pcs.	MP01820	■	■	■
G	Mechanical filter and HME TwinStar® HEPA, disposable, 55 mL (1.86 fl. oz.) deadspace, 50 pcs.	MP01801	■	■	■

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



BREATHING SYSTEM FILTERS AND HMEs

ELECTROSTATIC FILTERS

The CareStar® breathing system filter significantly contributes to infection prophylaxis in anesthesia and ventilation treatment through high retention rates.

- High-performance electrostatic filter
- Patient-side or device-side use
- Minimal workload and low cost
- Sampling connector with tethered cap
- Sampling port in convenient 45° angle
- Transparent housing for visual control
- Standardized connectors for safe connection to other components
- Clear labeling and red color coding for quick identification

A	Electrostatic filter CareStar® 45, disposable, 45 mL (1.52 fl. oz.) deadspace, 50 pcs.	MP01755	■	■
B	Electrostatic filter CareStar® 40A, disposable, 40 mL (1.35 fl. oz.) deadspace, angled, 50 pcs.	MP01765	■	■
C	Electrostatic filter CareStar® 30, disposable, 30 mL (1.01 fl. oz.) deadspace, 50 pcs.	MP01770	■	■

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A

MT-4387-2007



B

MT-4747-2007



C

MT-0011-2007



BREATHING SYSTEM FILTERS AND HMEs

MECHANICAL FILTERS

The SafeStar® mechanical high performance breathing system filter shows excellent retention rates, thereby significantly contributing to infection prophylaxis in anesthesia and ventilation treatment.

- High-performance mechanical HEPA filter
- High hydrophobicity
- Patient-side or device-side use
- Sampling connector with tethered cap
- Sampling port in convenient 45° angle
- Transparent housing for visual control
- Standardized connectors for safe connection to other components
- Clear labeling and red color coding for quick identification

A	Mechanical filter SafeStar® 80, disposable, 80 mL (2.71 fl. oz.) deadspace, 50 pcs.	MP01785	■	■
B	Mechanical filter SafeStar® 60A, disposable, 60 mL (2.03 fl. oz.) deadspace, angled, 50 pcs.	MP01795	■	■
C	Mechanical filter SafeStar® 55, disposable, 55 mL (1.86 fl. oz.) deadspace, 50 pcs.	MP01790	■	■

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



BREATHING SYSTEM FILTERS AND HMEs

HMEs

HumidStar® heat and moisture exchangers from Dräger partially take over the function of the upper airways and efficiently heat and humidify the inspired air.

- High return of heat and moisture
- Reduced condensation in the breathing circuit
- Minimal workload and low cost
- Sampling connector with tethered cap
- Sampling port in convenient 45° angle
- Transparent housing for visual control
- Standardized connectors for safe connection to other components
- Clear labeling and green color coding for quick identification

A	HME HumidStar® 55, disposable, 55 mL (1.86 fl. oz.) deadspace, 50 pcs.	MP01730	■	■	■
B	HME HumidStar® 25, disposable, 25 mL (0.85 fl. oz.) deadspace, 50 pcs.	MP01735	■	■	■
C	HME HumidStar® 10A, disposable, 10 mL (0.34 fl. oz.) deadspace, angled, 50 pcs.	MP01740	■	■	■
D	HME HumidStar® 2, disposable, 2 mL (0.07 fl. oz.) deadspace, 50 pcs.	MP01745	■	■	■

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



CATHETER MOUNTS AND AIRWAY CONNECTORS

4. CATHETER MOUNTS AND AIRWAY CONNECTORS

- Smoothbore products provide low resistance to flow while external coil reduces kinks and blockages
- Extendable product provides a lightweight and flexible solution
- Double swivel connector increases ease of use
- Double cap to maintain PEEP during suctioning or bronchoscopy
- Each catheter mount and airway connector is equipped with a red safety cap so that the system is completely closed and protected until directly prior to use on the patient
- 22F/15M connector on device side as a two-in-one solution
- Standardized connectors for safe connection to other components
- Luer Lock port for gas sampling
- Disposable convenience enhances workflow
- PVC-free design of gas leading components eliminates risk of softeners such as DEHP



A	Catheter mount ErgoStar® CM 40, disposable, 10 cm (3.94 in), double swivel elbow, 50 pcs.	MP01840	■	■	■
B	Catheter mount ErgoStar® CM 45, disposable, 15 cm (5.91 in), double swivel elbow, 50 pcs.	MP01845	■	■	■
C	Catheter mount ErgoStar® CM 50, disposable, 10 – 16 cm (3.94 – 6.3 in), extendable, 50 pcs.	MP01850	■	■	■
D	Catheter mount ErgoStar® CM 55, disposable, 15 cm (5.91 in), straight, 50 pcs.	MP01855	■	■	■
E	Catheter mount ErgoStar® CM 60, disposable, 10 cm (3.94 in), small, double swivel elbow, 50 pcs.	MP01860	■	■	■
F	Airway connector ErgoStar® AC 90, disposable, double swivel elbow w/ double cap, 100 pcs.	MP01890	■	■	■
G	Airway connector ErgoStar® AC 95, disposable, w/ Luer Lock sampling port, 100 pcs.	MP01895	■	■	■



COMBINATIONS OF FILTER AND TUBE EXTENTION

COMBINATIONS OF FILTER AND TUBE EXTENTION

A	CombiStar Filter HME flex	MP04240
	Combination of flexible tube extension MP01850 and filter TwinStar® 55 (MP01805).	
	CombiStar F-HME HEPA flex	MP04242
	Combination of flexible tube extension MP01850 and filter TwinStar® HEPA (MP01801).	
B	CombiStar Filter flex	MP04244
	Combination of flexible tube extension MP01850 and filter CareStar® 30 (MP01770).	
C	CombiStar mech. Filter flex	MP04246
	Combination of flexible tube extension MP01850 and filter SafeStar® 55 (MP01790).	
C	CombiStar A Filter HME	MP01671
	Combination of airway connector MP01845 and filter TwinStar® 55 (MP01805).	

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

COMBINATIONS OF FILTER AND TUBE EXTENTION

A	CombiStar Filter HME straight	MP04230
	Combination of straight tube extension MP01855 and filter TwinStar® 55 (MP01805).	
	CombiStar F-HME HEPA straight	MP04232
	Combination of straight tube extension MP01855 and filter TwinStar® HEPA (MP01801).	
B	CombiStar Filter straight	MP04234
	Combination of straight tube extension MP01855 and filter CareStar® 30 (MP01770).	
C	CombiStar mech. Filter straight	MP04236
	Combination of straight tube extension MP01855 and filter SafeStar® 55 (MP01790).	

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL


NIV MASKS AND ACCESSORIES

5. BREATHING MASKS

NIV MASKS AND ACCESSORIES

ClassicStar® NIV oronasal masks cover the nose and mouth of the patient. Therefore, they provide effective therapy even in case of mouth breathing. The mask design ensures high patient comfort combined with anatomical fit and effective seal due to individually adjustable mask cushion. The 6-point headgear provides a stable and secure fixation. Adjustable forehead support and forehead pad allow for distributing the pressure from the bridge of nose towards the forehead for more patient comfort. ClassicStar® NIV masks with standard elbow (SE) have been successfully tested for compatibility with Dräger ventilators which have NIV-option. Variants with Anit-Asphyxia Valve (AAV) are intended for Multi-Vendor-Accessory business. All ClassicStar® NIV masks are for single patient use. Sizing gauges for easy mask size identification as well as a pumpball are delivered together with the masks.

A	NIV oronasal mask ClassicStar® SE, disposable, size S	MP01573	■	■	■
	NIV oronasal mask ClassicStar® SE, disposable, size M	MP01574	■	■	■
	NIV oronasal mask ClassicStar® SE, disposable, size L	MP01575	■	■	■




NIV MASKS AND ACCESSORIES

Oronasal masks are often the first-line strategy in the initial management of acute respiratory failure as they cover nose and mouth. Nasal masks are often chosen if NIV therapy has to be prolonged due to improved comfort and reduced complications associated with oronasal masks. Similar to ClassicStar® oronasal masks, the nasal mask cushion can be individually adjusted to the patient's face providing high patient comfort and effective seal. The 4-point-headgear ensures a stable and secure fixation reducing mask displacement. All ClassicStar® NIV nasal masks are for single patient use. Sizing gauges for easy mask size identification as well as a pumpball are delivered together with the masks.

A	NIV nasal mask ClassicStar® non-vented, disposable, size S	MP01623	■	■
	NIV nasal mask ClassicStar® non-vented, disposable, size M	MP01624	■	■
	NIV nasal mask ClassicStar® non-vented, disposable, size L	MP01625	■	■

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



NIV MASKS AND ACCESSORIES

The ClassicStar® NIV total-face mask seals around the perimeter of the face, covering mouth, nose and eyes. These areas are less pressure sensitive and have smoother facial contours compared to the areas a oronasal masks covers (nasal bridge). The mask is suitable for patients with nasal bridge breakdown or facial irregularities. The air filled cushion can be adjusted individually to the patients face, providing high patient comfort and effective seal.

The 5-point headgear which includes a headstrap at the top of the mask provides a stable and secure attachment, thus preventing mask displacement and reducing pressure. Quick release clips are integrated for quick removal of the mask.



EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

**A NIV total-face mask ClassicStar® SE,
disposable, size S**

MP02910

**NIV total-face mask ClassicStar® SE,
disposable, size M**

MP02911

**NIV total-face mask ClassicStar® SE,
disposable, size L**



MP02912

NIV MASKS AND ACCESSORIES

The ClassicStar® plus mask series features a soft and anatomically shaped silicone lip as sealing interface to the patients face. The mask has an easy to adjust forehead pad with silicone cushion which allows optimal adaption to the sensitive nasal bridge. The mask is available in four sizes (S, M, L, XL) as oronasal and nasal mask. ClassicStar® plus NIV masks with standard elbow (SE) have been successfully tested for compatibility with Dräger ventilators which have NIV-option. Variants with Anti-Asphyxia-Valve (AAV) are intended for Multi-Vendor-Accessory business.

A	NIV oronasal mask ClassicStar® plus SE, disposable, size S	MP04711			
	NIV oronasal mask ClassicStar® plus SE, disposable, size M	MP04712			
	NIV oronasal mask ClassicStar® plus SE, disposable, size L	MP04713			
	NIV oronasal mask ClassicStar® plus SE, disposable, size XL	MP04714			
B	NIV nasal mask ClassicStar® plus non vented, disposable, size S	MP04741			
	NIV nasal mask ClassicStar® plus non vented, disposable, size M	MP04742			
	NIV nasal mask ClassicStar® plus non vented, disposable, size L	MP04743			
	NIV nasal mask ClassicStar® plus non vented, disposable, size XL	MP04744			

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



NIV MASKS AND ACCESSORIES

NovaStar® TS NIV oronasal masks cover the nose and mouth of the patient. Therefore, they provide effective therapy even in case of mouth breathing.

The NovaStar® TS NIV masks have very soft silicone-gel filled mask cushion. In combination with the pliable ring embedded inside the flexible, clear mask body, – which allows the mask to be bent and adjusted individually to the patient's face ('customized fit'), – the masks provide ultimate comfort and effective seal. The forehead support can be adapted gradually to the patient's forehead providing stable mask position and enhanced comfort. The unique magnetic clips support fast and easy fixation.

All NovaStar® TS NIV masks are for multi patient use. They can be disinfected thermally with hot steam or chemically (manually). Replacement headgears and clips are available individually (cf. accessories for NIV masks).



EVITA® INFINITY® V500
EVITA® V300
EVITA® XL




A NIV oronasal mask NovaStar® TS SE, reusable, size S	<u>MP01579</u>	■ ■ ■
NIV oronasal mask NovaStar® TS SE, reusable, size M	<u>MP01580</u>	■ ■ ■
NIV oronasal mask NovaStar® TS SE, reusable, size L	<u>MP01581</u>	■ ■ ■

NIV MASKS AND ACCESSORIES

ACCESSORIES FOR NIV MASKS

A	Headgear for NovaStar® TS SE NIV oronasal mask	MP01583
B	Magnetic clip for headgears for NovaStar® TS SE and AAV, 4 pcs.	MP01585
C	Sizing gauge for ClassicStar® SE/AAV NIV full-face masks, 3 pcs.	MP01592
D	Sizing gauge for ClassicStar® non-vented/vented NIV nasal masks, 3 pcs.	MP01632
E	Two-sided pump ball for air management	MP01590


EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



ANESTHESIA MASKS AND ACCESSORIES

ANESTHESIA MASKS AND ACCESSORIES

A	① Anesthesia mask LiteStar®, disposable, w/o hook ring, size 1, prenat., 30 pcs.	MP01501			
	② Anesthesia mask LiteStar®, disposable, w/o hook ring, size 2, baby, 30 pcs.	MP01502			
	③ Anesthesia mask LiteStar®, disposable, w/o hook ring, size 3, infant, 30 pcs.	MP01503			
	④ Anesthesia mask LiteStar®, disposable, w/o hook ring, size 4, adult S, 30 pcs.	MP01504			
	⑤ Anesthesia mask LiteStar®, disposable, w/o hook ring, size 5, adult M, 30 pcs.	MP01505			
	⑥ Anesthesia mask LiteStar®, disposable, w/o hook ring, size 6, adult L, 30 pcs.	MP01506			
B	① Hook ring for LiteStar® anesthesia masks, disposable, 15 mm (0.59 in) OD, for sizes 1 – 2, 30 pcs.	MP01551			
	② Hook ring for LiteStar® anesthesia masks, disposable, 22 mm (0.87 in) ID, for sizes 3 – 6, 30 pcs.	MP01552			

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

ANESTHESIA MASKS AND ACCESSORIES

A	① Anesthesia mask ComfortStar®, disposable, hook rg, std., size 0, prenat., 20 pcs.	MP01510
	② Anesthesia mask ComfortStar®, disposable, hook rg, std., size 1, neo., 20 pcs.	MP01511
	③ Anesthesia mask ComfortStar®, disposable, hook rg, std., size 2, baby, 20 pcs.	MP01512
	④ Anesthesia mask ComfortStar®, disposable, hook rg, std., size 3, inf., 20 pcs.	MP01513
	⑤ Anesthesia mask ComfortStar®, disposable, hook rg, std., size 4, adt. S, 20 pcs.	MP01514
	⑥ Anesthesia mask ComfortStar®, disposable, hook rg, std., size 5, adt. M, 20 pcs.	MP01515
	⑦ Anesthesia mask ComfortStar®, disposable, hook rg, std., size 6, adt. L, 20 pcs.	MP01516






EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

ANESTHESIA MASKS AND ACCESSORIES

ACCESSORIES FOR ANESTHESIA MASKS

A	① Headgear, reusable, silicone, small	MP01557
	② Headgear, reusable, silicone, large	MP01558
B	① Headgear, disposable, soft fabric, small	MP01559
	② Headgear, disposable, soft fabric, large	MP01560
C	① Headgear, disposable, soft rubber, small	MP01561
	② Headgear, disposable, soft rubber, large	MP01562

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



PACK2GO®/SET2GO

6. PACK2GO®

Dräger has developed a comprehensive range of disposable accessories. These accessories can now also be ordered together in a Pack2Go® package, in order to benefit the hospital's costs and logistical/clinical processes. Pack2Go® package solutions include dedicated Dräger and general accessories, which can be used as a system solution with Dräger devices.

A	Pack2Go® Evita® invasive vent. pack w/ breathing circuit VentStar®, basic, + cuvette	MP01656
B	Pack2Go® Evita® invasive vent. pack w/ breathing circuit VentStar®, basic, 1.8 m	MP01657
C	Pack2Go® Evita® invasive vent. pack w/ breathing circuit VentStar®, coaxial, + cuvette	MP01676
D	Pack2Go® Evita® invasive vent. pack w/ breathing circuit VentStar®, coaxial, 1.8 m	MP01677

SET2GO

With Dräger's Set2go solutions we offer you totally customised sets consisting of an individual selection of single-patient-use accessories. Dräger Set2go makes using and ordering consumables and accessories much simpler for you. You just need one article number to get all the items in a single package and accomplish your specific task.

Please contact your local Dräger sales organisation to create your individual accessory solution.

EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



RESUSCITATORS

7. RESUSCITATORS

The MR-100 is a reusable, silicone, manual resuscitator that is so compact it can be folded to approximately half its normal size when not in use. A pop-off valve is included, to prevent high airway pressures. The resuscitator is available in three sizes: adult (A), child (P) and infant (N).

A	Man. breathing bag MR-100, reusable, latex-free, A, w/ pressure valve, 60±10 mbar	5702881
B	Man. breathing bag MR-100, reusable, latex-free, P, w/ pressure valve, relief >40±5 mbar	5702321
C	Man. breathing bag MR-100, reusable, latex-free, N, w/ pressure valve, relief >40±5 mbar	5702322

ACCESSORIES

D	O ₂ -Reservoir MR-100, 2,500 mL*	5703701
E	O ₂ -Reservoir MR-100, 600 mL*	5703702

*Made by GaleMed Corporation, Taiwan



Not all products, features, or services are for sale in all countries.

Mentioned Trademarks are only registered in certain countries and not necessarily in the country in which this material is released. Go to www.draeger.com/trademarks to find the current status.

CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany

www.draeger.com

Manufacturer:

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany

REGION DACH

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION EUROPE

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION MIDDLE EAST, AFRICA

Drägerwerk AG & Co. KGaA
Branch Office
P.O. Box 505108
Dubai, United Arab Emirates
Tel +971 4 4294 600
Fax +971 4 4294 699
contactuae@draeger.com

REGION ASIA / PACIFIC

Draeger Singapore Pte. Ltd.
25 International Business Park
#04-20/21 German Centre
Singapore 609916
Tel +65 6308 9400
Fax +65 6308 9401
asia.pacific@draeger.com

REGION NORTH AMERICA

Draeger, Inc.
3135 Quarry Road
Telford, PA 18969-1042, USA
Tel +1 800 4DRAGER
(+1 800 437 2437)
Fax +1 215 723 5935
info.usa@draeger.com

REGION CENTRAL AND SOUTH AMERICA

Dräger Panama S. de R.L.
59 East Street, Nuevo Paitilla,
House 30, San Francisco Town
Panama City, Panama
Tel +507 377 9100
Fax +507 377 9130
servicioalcliente@draeger.com

Locate your Regional Sales
Representative at:
www.draeger.com/contact



