



**ResMed**

**Mobi™**

Oxygen is life  
*Mobi is my reassurance*



[ResMed.com/GoMobi](https://www.resmed.com/GoMobi)



# Get the most out of life

ResMed Mobi is a portable oxygen concentrator for everyday living. It features a lightweight design, long battery life and delivers the oxygen you need to live your life.



# Go further

Mobi is designed to help you enjoy an active lifestyle. For long-distance travel, Mobi can go where you go.

## Oxygen for you

Mobi uses Pulse-Wave™ technology to efficiently deliver your selected dose of oxygen each time you inhale. As your activity level changes, Mobi senses your breathing and automatically switches between active and rest mode.

## Simple to use

Mobi is designed for one-touch startup and simple setting changes. An easy-to-read display shows battery life and indicates when each dose of oxygen is delivered.



\*Accessories sold separately.

Active lifestyle. From daily errands to  
where you go.

## Casual style

The Mobi carry bag is both stylish and functional. It features a casual but sophisticated design that conceals tubing and batteries. Or choose the Mobi backpack – great for day trips and travel.

## A name you can trust

ResMed has been changing lives with every breath for nearly thirty years. Mobi is the latest addition to our award-winning family of respiratory solutions.



# Mobi™

Portable oxygen concentrator



# ResMed

## Mobi accessories



### Mobi external battery

Offering up to 4.5 hours of additional running time at setting 1, the Mobi external battery allows you to get out and enjoy life longer.



### Mobi backpack

The Mobi backpack keeps your machine, batteries and nasal cannula secure and out of sight, while the backpack straps conveniently conceal the cannula.



### Mobi mounting straps (4)

A set of four adjustable straps feature hook and clip attachments to help you securely attach Mobi to a scooter, walker or wheelchair.



### Nasal cannula

7-ft silicone nasal cannula offers a comfortable fit and ample length for freedom of movement.

\*Accessories sold separately.

## Technical specifications

### Mobi device

Weight	5.5 lbs (2.49 kg)		
Size (W x H x D)	9.60" x 7.95" x 4.49" (244 mm x 202 mm x 114 mm)		
<b>Oxygen delivery</b>	Pulse oxygen therapy		
Device settings	1, 2, 3, 4		
<b>Battery life*</b>	<b>Setting</b>	<b>Internal battery</b>	<b>External</b>
	1	8.6 hrs	4.5 hrs
	2	6 hrs	3 hrs
	3	3.8 hrs	1.9 hrs
	4	3.2 hrs	1.6 hrs

\*Hours are approximate and based on a breath rate of 20 bpm. The device will operate up to the battery life (hours) specified.

## ResMed.com/Mobi

Distributed by ResMed Corp, 9001 Spectrum Center Boulevard, San Diego, CA 92123 USA +1 858 836 5000 or 1 800 424 0737 (toll free). See ResMed.com for other ResMed locations worldwide. Mobi and Pulse-Wave are trademarks and/or registered trademarks of the ResMed family of companies. Specifications may change without notice. For patent and other intellectual property information, see ResMed.com/ip. © 2018 ResMed. 10110452/2 2018-05